aquaticPRO

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DESCRIPTION

A waterproof protective cover designed to keep casts, elastic bandages, wounds, and medical dressings dry while showering or bathing. Provides a watertight seal without the need for additional adhesives bandages.

INSTRUCTIONS FOR USE

1. Gently stretch the silicone opening to create enough space for inserting the limb.
2. Carefully insert the limb with the cast or bandage into the protector.
3. Ensure that the cast or bandage is completely covered by the protector.
4. Check that the silicone seal fits properly against the skin without gaps or folds, and smooth out any wrinkles in the film for better waterproof sealing.
5. Carefully remove the protector by gently pulling the silicone opening. Do not use excessive force.
6. Wipe and allow the product to air dry naturally after use.

SAFETY WARNINGS

- Do not allow direct contact between the protector and open wounds.

- Use only for temporary protection during exposure to water.
- Do not use the product for swimming or diving.
- Exercise caution when walking on wet or slippery surfaces, as the protector may become slippery.
- Keep the product away from fire or high temperatures.
- Adult supervision is recommended when used by children. There is a potential risk of choking or suffocation.
- Avoid contact with sharp objects, nails, scissors, or hard surfaces that may damage or tear the silicone seal and the protective cover.
- Do not use the product if it shows cracks, holes, or signs of damage.

CLEANING & MAINTENANCE

Clean with lukewarm water and mild soap. Do not use chemical cleaners or bleach. Allow the product to air dry naturally before storage. Store in a cool and dry place away from direct sunlight.

IMPORTANT INFORMATION

The product is reusable. For external use only. If irritation or discomfort occurs, discontinue use and consult a physician.

TECHNICAL SPECIFICATIONS

Materials: PVC, Silicone - Dimensions: 350 x 240 mm

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